

FIL SSF

lente a fissazione sclerale

scleral fixation lens

sutureless



REVIEW

PREMIUM IOL

FIL SSF

Sutureless Scleral Fixation


Lente dal design innovativo che, con i caratteristici arpioni, garantisce l'ancoraggio al tessuto sclerale senza l'applicazione di punti di sutura, ne derivano semplicità operativa e tempi di intervento drasticamente ridotti.

Lens with innovative design that assures the anchoring to the scleral tissue by means of special harpoons without the need of suture stitches, leading to a simpler operation and drastically reduced surgery time.

FIL SSF

MONOFOCAL

Diametro ottico Optic diameter	6.5 mm
Diametro totale Total diameter	13.2 mm
Angolazione ansa Haptic angulation	10°
Materiale Material	foldable acrylic with 25% H ₂ O and UV filter
Indice di rifrazione Refractive index	1,461 (546 nm , 20°C in water)
Gamma poteri Diopter range	from -5.00 to +35.00 (step 0.5D)
Costante A consigliata Recommended A constant	118.5
Iniettore consigliato Recommended injector system	Medicel Viscojet 2.2 or 2.7




*Posizionamento corretto
*Correct positioning

FIL SSF

TORIC

Diottrie disponibili* (sfera) Available diopters* (sphere)	Cilindri disponibili* (cil) Available cylinders* (cyl)
step 0.5 D	step 0.25 D axis 0-180, step 1°
from -5.00 to +9.50 D	from 1.00 to 10.00 D
from +10.00 to +19.50 D	from 1.00 to 15.00 D
from +20.00 to +25.00 D	from 1.00 to 10.00 D
from +25.50 to +35.00 D	from 1.00 to 6.00 D
Costante A consigliata Recommended A constant	118.5
Iniettore consigliato Recommended injector system	Medicel Viscojet 2.2 or 2.7

Per il calcolo <http://www.soleko-iol.it/iol-calculator/>



*Posizionamento corretto
*Correct positioning

Real Axis Technology
L'asse del cilindro viene impostato in laboratorio al momento della costruzione della lente ed è indicato sul disegno tecnico allegato. La lente deve essere sempre posizionata a 0°-180°.
The cylinder axis is set during the construction and it is reported on the technical drawing belonging the lens. The lens must always be positioned at 0°-180°.

*La lente viene etichettata riportato l'equivalente sferico (SE = sfera + 1/2 cil). Per sfere o cilindri fuori gamma contattare ioldivision@soleko.it
*On the lens label it is reported the spherical equivalent (SE = sph + 1/2 cyl). For out of range spheres or cylinders please contact ioldivision@soleko.it